

A Case Study Utilizing Caris GPSai To Ensure Correct Diagnosis Of Primary Tumor Site

Background:

- Caris MI Profile[™] Tissue testing performed at Caris Life Sciences[®] includes Whole Exome Sequencing (WES), Whole Transcriptome Sequencing (WTS) and a panel of Immunohistochemical stains (IHCs). In addition, Caris GPSai[™] scores are generated for cases that are successful sequenced for WES and WTS.
- **GPSai is an Artificial Intelligence (AI) algorithm** that compares the WES and WTS data generated to our library of 294,000+ molecular database of comprehensive tumor profile to generate a score predicting the primary tumor site.

• Clinical utility of GPSai includes:

- Identify the tissue of origin for cancer of unknown primary (CUP) patients.
- Identify patients that are incorrectly diagnosed.

Presented Case

• Clinical presentation:

- 50 Year old female.
- Presented with bilateral adnexal masses (largest measuring 22cm in greatest dimension).
- Underwent total abdominal hysterectomy, bilateral salpingo-oophorectomy with pelvic debulking of tumor.

• The pathological diagnosis:

- Ovarian endometrioid adenocarcinoma FIGO grade 1.
- This was based on morphological examination, no IHCs done.
- Specimen received and testing performed at Caris.
 - A representative specimen of the right adnexal lesion was sent to Caris to help decide the best options for therapy.
 - MI Profile Tissue for ovarian tumor was run, which includes WES and WTS and a panel of IHCs including ER, PR and MMR (MLH1, MSH2, MSH6 and PMS2).
 - Based on the sequencing results for WES and WTS, GPSai score was generated. This highly favored colorectal cancer (CRC).
 - Subsequently, a panel of diagnostic IHCs were performed that showed the tumor to be positive for CK20, CDX2 and SATB-2 and negative for CK7, PAX-8 and ER. This immunohistochemical profile also strongly favored CRC.
- The findings were discussed in detail with ordering oncologist.

Key Findings: Specimens



(A) Moderately differentiated adenocarcinoma (H&E 20X). Positive staining of CK20 (B), CDX2 (C) and SATB2 (D). Negative staining of CK7 (E), PAX8(F), ER (G) and PR (H).

- Patient underwent colonoscopy that showed large nearly obstructive colonic lesion.
- This case was reviewed at the respective hospital tumor board and there was consensus that this patient is best managed as metastatic colonic adenocarcinoma.